



5

## **IHE Technical Frameworks General Introduction**

10

### **Appendix A: IHE Actor Definitions**

15

20

**Revision 1.0  
July 1, 2014**

25

**Please verify you have the most recent version of this document, which is published [here](#).**

## IHE Actor Definitions

Actors are information systems or components of information systems that produce, manage, or act on information associated with operational activities in the enterprise.

30

<b>Archive</b>	Stores the radiotherapy (RT) DICOM SOP Classes in addition to the computed tomography (CT) images and is capable of transmitting them.
<b>Acquisition Modality</b>	Acquires medical images, waveforms or measurements while a patient or specimen is present (e.g., a Computed Tomography scanner, a specimen microscope or a hemodynamic measurement system).
<b>Acquisition Modality Importer</b>	Integrates a non-DICOM-ready modality into workflows.
<b>ADT/Patient Registration</b>	Adds or updates patient demographic and encounter information, (Admission/Discharge/Transfer) to a Practice Management System (PMS) or a Hospital Information System (HIS), for example.
<b>Analyzer</b>	Performs testing on biological specimens and reports the resulting observations.
<b>Audit Repository</b>	Stores audit events.
<b>Automation Manager</b>	Schedules and manages operations on biological specimens on a set of robotic devices in a clinical laboratory.
<b>Charge Processor</b>	Receives posted charges and serves as a component of the financial system (for instance a PMS or billing system).
<b>Client Authentication Agent</b>	Provides local management of user authentication.
<b>Code Set Master</b>	Manages and provides code sets for laboratory panels, tests and observations
<b>Code Set Consumer</b>	Maintains its internal dictionary of laboratory panels, tests and observations relying on an external code set provided by a Code Set Master.
<b>Content Consumer</b>	Views, imports, or processes content created by a Content Creator Actor.
<b>Content Creator</b>	Creates content for later consumption by a Content Consumer.
<b>Context Manager</b>	Brokers context information (such as current patient subject) between context.
<b>Contourer</b>	Consumes CTs and creates RT Structure Sets.

<b>Department System Scheduler/Order Filler</b>	Manages the scheduling and performance of orders for diagnostic services (e.g., cardiology, radiology, cytology & anatomic pathology examinations, in vitro diagnostic testing) in a department or setting offering such services.
<b>Device Observation Consumer</b>	Consumes Patient Care Device data.
<b>Device Observation Filter</b>	Provides publish/subscribe-based filtering of Patient Care Device data.
<b>Device Observation Reporter</b>	Receives data from patient care devices and makes it available with consistent syntax and semantics.
<b>Display</b>	Requests and displays preformatted (“presentation-ready”) data.
<b>DNS Server</b>	Provides authoritative network location information.
<b>Document Consumer</b>	Queries for document metadata meeting certain criteria, and may retrieve selected documents
<b>Document Recipient</b>	Receives a set of documents (e.g., a viewer or Health Record system).
<b>Document Registry</b>	Manages and provides query access to meta-data about registered documents such as the document type, subject and where it is stored.
<b>Document Repository</b>	Persistently stores documents and makes them available.
<b>Document Source</b>	Produces and publishes documents.
<b>Dose Displayer</b>	Consumes a Dosimetric Plan, CT, Structure Set and an RT Dose and displays the dose.
<b>Dosimeter Planner</b>	Consumes CT, and RT Structure Set, a Geometric Plan, and creates a Dosimetric Plan and an RT Dose.
<b>Enterprise Report Repository</b>	Receives Structured Report Export Transactions from the Report Manager and stores the report.
<b>Evidence Creator</b>	Creates evidence objects such as derived images, presentation states, Key Image Notes, and/or measurements (Evidence Documents), and transmits them to an Image Archive.
<b>Export Manager</b>	De-identifies, pseudonymizes and exports imaging data.
<b>Export Selector</b>	Selects and tags imaging data for export, e.g. to a teaching file or a clinical trial.
<b>External Report Repository Access</b>	Retrieves clinical reports generated outside the imaging department.

<b>Geometric Planner</b>	Consumes a CT and RT Structure Set and creates a Geometric Plan.
<b>Image Archive</b>	Provides long term storage of imaging data such as images, measurements, presentation states, and manifests (e.g., a PACS).
<b>Image Display</b>	Queries for, retrieves and displays imaging data.
<b>Image Manager</b>	Manages and provides access to stored imaging objects.
<b>Imaging Document Consumer</b>	Retrieves imaging data based on references in a retrieved imaging manifest.
<b>Imaging Document Source</b>	Publishes imaging data and makes it available for retrieval.
<b>Implantable Device Cardiac Consumer</b>	Receives cardiac implantable device data.
<b>Implantable Device Cardiac Reporter</b>	Reports data from Cardiac Implantable Devices.
<b>Importer</b>	Imports and reconciles imaging data.
<b>Information Source</b>	Provides presentation-ready information or documents.
<b>Initiating Gateway</b>	Supports all outgoing inter-community communications.
<b>Infusion Order Consumer</b>	Programs an infusion pump based on order information.
<b>Infusion Order Programmer</b>	Provides infusion pump order information.
<b>Kerberos Authentication Server</b>	Provides central authentication of enterprise users.
<b>Kerberized Server</b>	Receives user authentication information
<b>Label Broker</b>	Delivers pre-identified labels or labeled containers for the collection of biological specimens requested for the performance of in vitro diagnostic testing orders.
<b>Label Information Provider</b>	Manages and provides specimen labeling instructions.
<b>Master Patient Index (MPI)</b>	Maintains unique enterprise-wide identifiers for patients.
<b>Order Filler</b>	See Department System Scheduler/Order Filler
<b>Order Placer</b>	Generates and distributes orders to various departments.
<b>Order Result Tracker</b>	Stores and manages laboratory or pathology observations and associated results.

<b>Patient Context Participant</b>	Initiates and responds to patient context changes.
<b>Patient Demographics Consumer</b>	Uses demographic information.
<b>Patient Demographics Supplier</b>	Provides a searchable repository of patient demographic information.
<b>Patient Encounter Consumer</b>	Uses patient encounter information.
<b>Patient Encounter Supplier</b>	Manages and provides information about new and existing encounters.
<b>Patient Identifier Cross-reference Consumer</b>	Obtains and uses the identification of a patient in different Patient Identifier Domains.
<b>Patient Identifier Cross-reference Manager</b>	Provides cross-referencing of patient identifiers across Patient Identifier Domains.
<b>Patient Identity Source</b>	Provides a unique identifier and identity traits for each patient.
<b>Performed Procedure Step Manager</b>	Re-distributes Performed Procedure Step information.
<b>Personnel White Pages Consumer</b>	Uses information from the Personnel White Pages Directory.
<b>Personnel White Pages Directory</b>	Holds authoritative Personnel White Pages information.
<b>Pre/Post-Processor</b>	Processes biological specimens (e.g., sorters, aliquoters, decappers, specimen conveyors, specimen storage systems).
<b>Point Of Care Data Manager</b>	Manages a set of point of care in vitro diagnostic devices and centralizes their results.
<b>Point Of Care Results Generator</b>	Produces in vitro diagnostic observations on the point of care.
<b>Portable Media Creator</b>	Marshals content and writes it to physical media such as CDs and USB storage.
<b>Portable Media Importer</b>	Reads imaging data on media, and reconciles and stores the instances.
<b>Post-Processing Manager</b>	Schedules and manages post-processing worklists.
<b>Print Composer</b>	Generates DICOM print requests to a Print Server.
<b>Print Server</b>	Renders hardcopy images based on print requests.
<b>Receiver</b>	Receives exported instances over the network.

---

<b>Report Creator</b>	Generates and transmits diagnostic reports.
<b>Report Manager</b>	Manages reporting worklists and short-term storage of Report objects during the reporting process then distributes reports to report repositories.
<b>Report Reader</b>	Accesses reports through network query/retrieve or reading interchange media and allow the user to view reports presented as DICOM Structured Reporting Objects.
<b>Report Receiver</b>	Receives and stores anatomic pathology reports.
<b>Report Repository</b>	Provides long-term storage and retrieval of diagnostic reports.
<b>Report Sender</b>	Creates or has access to anatomic pathology reports in final status and sends these reports or parts of these reports.
<b>Requester</b>	Forwards a part of a laboratory testing requisition together with the associated biological specimens to a subcontracting performer, and collects the resulting observations.
<b>Responding Gateway</b>	Supports all incoming inter-community communications.
<b>Secure Application</b>	Complies with IHE rules for user authentication, secure communications, security audit recording, and security policy enforcement. The specified application is compliant.
<b>Secure Node</b>	Complies with IHE rules for user authentication, secure communications, security audit recording and security policy enforcement. All software on the system is compliant.
<b>Subcontractor</b>	Manages the scheduling and performance of in vitro diagnostic tests in a subcontracting laboratory, and reports the resulting observations to the requesting laboratory.
<b>Time Client</b>	Maintains system clock synchronization with UTC based on time synchronization with one or more Time Servers.
<b>Time Server</b>	Provides NTP time services to Time Clients.
<b>User Context Participant</b>	Initiates and responds to user context changes.
<b>X-Service Provider</b>	Provides a service that needs an X-User Assertion.
<b>X-Service User</b>	Makes service requests of an X-Service Provider.